



# Aspect™ Systems Installation Quick Start

**Aspect Systems Installation Quick Start**

**© INFOTECH Soft, Inc.**

9200 South Dadeland Boulevard • Suite 620  
Miami • Florida • 33156

Phone 305.670.5111 • Fax 305.670.7722

Email [solutions@infotechsoft.com](mailto:solutions@infotechsoft.com)

Web <http://www.infotechsoft.com>

Last Updated September 30, 2014

The information of this document may change without notice.

In case of any problem with this document, please send an email to [solutions@infotechsoft.com](mailto:solutions@infotechsoft.com).

Aspect™, Aspect Rating Station, Aspect Control Center, Aspect Data Source, and ClinicalSync™ are trademarks of INFOTECH Soft, Inc.

## Aspect Systems Hardware Setup

The Aspect System can be configured as an independent network, as part of a campus network, or as systems distributed across the Internet. In all cases, **the Control Center computer must have a static IP** in order for the Rating Stations to connect and synchronize data.

### Sample Hardware for the Aspect System


---

The following chart lists examples of available computers and network equipment that meet the minimum requirements of the Aspect System. For more information about the specific requirements of the Aspect Control Center and Rating Station, see the *Aspect Systems Installation Guide*.

<u>Function</u>	<u>Minimum Requirements</u>	<u>Example</u>
Control Center	512MB RAM, 30 GB hard drive, 10/100 Ethernet	Desktop: Dell Dimension 2400
Rating Station	512MB RAM, 10 GB hard drive, 10/100 Ethernet	Tablet: HP tc1100, Motion M1300 Laptop: Dell Latitude 100L Desktop: Dell Dimension 2400
Network Switch	4-Port 10 Mbps	Netgear 10/100 4-port switch, Gigafast 10/100 8-port switch, 3Com 10/100 4-port switch

## Aspect Control Center Installation

The following flow chart illustrates the installation process for the Aspect Control Center. Detailed instructions for installing the Aspect Control Center can be found in the *Aspect Systems Installation Guide*.



## Aspect Rating Station Installation

The following flow chart illustrates the installation process for the Aspect Rating Station. Detailed instructions for installing the Aspect Rating Station can be found in the *Aspect Systems Installation Guide*.

